

Product Datasheet

NTNG2 antibody (orb67221)

Description	Rabbit polyclonal to NTNG2.Netrin-G1 and netrin-G2, belonging to a vertebrate-specific subfamily of ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human NTNG2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	60
Uniprot ID	Q96CW9
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.

