


Product Datasheet

TDG antibody (orb67223)

Description	Rabbit polyclonal to TDG.G/T mismatch-specific thymine DNA glycosylase (TDG) belongs to the TDG/mug ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human TDG
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	46
Uniprot ID	Q13569
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...