


Product Datasheet

Histone H2B antibody (orb672294)

Description	Histone H2B antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from the Internal region of Human Histone H2B.
Target	HIST1H2BB/HIST1H2BC/HIST1H2BD/HIST1H2BH/HIST1H2BJ/HIST1H2BK/HIST1H2BL/HIST1H2BM/HIST1H2BN/HIST1H2BO
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P58876, P62807, P33778, Q99877, P23527, P06899, O60814, Q93079, Q99879, Q99880
Dilution Range	WB: 1:500-1:2000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using Histone H2B antibody (orb672294).