



Product Datasheet

TET1 antibody (orb67230)

Description	Rabbit polyclonal to TET1. TET1 is a methylcytosine dioxygenase that catalyzes the conversion of the ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, WB
Immunogen	C-terminal region of human TET1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	235
Uniprot ID	Q8NFU7
Dilution Range	WB:1:500-1:1000, ChIP: 1:100- 1:200
Expiration Date	12 months from date of receipt.



Chromatography analysis of human Jurkat ...



Western blot analysis of total cell extr...