


Product Datasheet

TET1 antibody (orb67232)

Description	Rabbit polyclonal to TET1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of mouse TET1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	235
Uniprot ID	Q8NFU7
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...