



Product Datasheet

TET1 antibody pair (orb67233)

Description	Rabbit polyclonal to TET1. 50 μl TET1 (N-term1) &50 μl TET1 (C-term2) .vPair antibodies repres...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	-
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	235
Uniprot ID	Q8NFU7
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...