


Product Datasheet

HP1BP3 antibody (orb67240)

Description	Rabbit polyclonal to HP1BP3. Heterochromatin protein 1-binding protein 3 (HP1BP3) is a component of h...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human HP1BP3
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	61
Uniprot ID	Q5SSJ5
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...