


Product Datasheet

TACR1 antibody (orb67248)

Description	Rabbit polyclonal to TACR1.Substance-P receptor /TACR1 is a multi-pass membrane protein that belongs...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human TACR1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	46
Uniprot ID	P25103
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...