


Product Datasheet

ZNF185 antibody (orb67251)

Description	Rabbit polyclonal to ZNF185. Zinc finger protein 185 (ZNF185) may function as a tumor-suppressor prot...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human ZNF185
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	74
Uniprot ID	O15231
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...