


Product Datasheet

TRIM33 antibody (orb67253)

Description	Rabbit polyclonal to TRIM33. Tripartite motif-containing protein 33 (TRIM33) is a nuclear protein loc...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human TRIM33
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	123
Uniprot ID	Q9UPN9
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...