


Product Datasheet

SHC3 antibody (orb67257)

Description	Rabbit polyclonal to SHC3.SHC-transforming protein 3 (SHC3) is a signaling adapter that couples acti...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SHC3
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	64
Uniprot ID	Q92529
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...