


Product Datasheet

STAT5B antibody (orb67261)

Description	Rabbit polyclonal to STAT5B. Signal transducer and activator of transcription 5B (STAT5B) belongs to ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human STAT5B
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	90
Uniprot ID	P51692
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...