


Product Datasheet

PAK2 antibody (orb67263)

Description	Rabbit polyclonal to PAK2.p21-activated kinase 2 (PAK2) belongs to the protein kinase superfamily, S...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human PAK2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	58
Uniprot ID	Q13177
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...