


Product Datasheet

SOCS6 antibody (orb67266)

Description	Rabbit polyclonal to SOCS6. Suppressor of cytokine signaling 6 (SOCS6) is a member of the SOCS family...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SOCS6
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	60
Uniprot ID	O14544
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...