


Product Datasheet

MMP20 antibody (orb67268)

Description	Rabbit polyclonal to MMP20. MMP20 is a matrix metalloproteinase belonging to the peptidase M10A famil...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human MMP20
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	54
Uniprot ID	O60882
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...