


Product Datasheet

TLK1 antibody (orb67284)

Description	Rabbit polyclonal to TLK1. Serine/threonine-protein kinase tousel-like 1 (TLK1) belongs to the prote...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human TLK1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	87
Uniprot ID	Q9UKI8
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...