


Product Datasheet

NRIP1 antibody (orb67285)

Description	Rabbit polyclonal to NRIP1. Nuclear receptor-interacting protein 1 (NRIP1) is also known as receptor-...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human NRIP1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	127
Uniprot ID	P48552
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...