


## Product Datasheet

### PIK3AP1 antibody (orb67286)

|                            |                                                                                                          |
|----------------------------|----------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal to PIK3AP1. Phosphoinositide 3-kinase adapter protein 1 (PIK3AP1) is involved in th... |
| <b>Species/Host</b>        | Rabbit                                                                                                   |
| <b>Reactivity</b>          | Human, Mouse                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                             |
| <b>Tested Applications</b> | IP, WB                                                                                                   |
| <b>Immunogen</b>           | N-terminal region of human PIK3AP1                                                                       |
| <b>Preservatives</b>       | 0.01% NaN <sub>3</sub>                                                                                   |
| <b>Storage</b>             | -20°C                                                                                                    |
| <b>Note</b>                | For research use only                                                                                    |
| <b>Clonality</b>           | Polyclonal                                                                                               |
| <b>Purity</b>              | Affinity purification                                                                                    |
| <b>MW</b>                  | 90                                                                                                       |
| <b>Uniprot ID</b>          | <b>Q6ZUJ8</b>                                                                                            |
| <b>Dilution Range</b>      | WB:1:1000-1:3000                                                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                          |



Western blot analysis of total cell extr...