


Product Datasheet

PRKCA antibody (orb67290)

Description	Rabbit polyclonal to PRKCA. Protein kinase C alpha type (PRKCA) is a peripheral membrane protein belo...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human PRKCA
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	77
Uniprot ID	P17252
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...