

Product Datasheet

MX2 antibody (orb672902)

Description	MX2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Synthesized peptide derived from the Internal region of Human MxB.
Target	MX2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P20592
Dilution Range	ELISA: 1:20000
Expiration Date	12 months from date of receipt.