


Product Datasheet

PTPN6 antibody (orb67294)

Description	Rabbit polyclonal to PTPN6. Tyrosine-protein phosphatase non-receptor type 6 (PTPN6) is also known as...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human PTPN6
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	68
Uniprot ID	P29350
Dilution Range	WB:1:5000-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...