


Product Datasheet

NKX2 antibody (orb67297)

Description	Rabbit polyclonal to NKX2.Homeobox protein Nkx-2.3 (NKX2-3) is a member of the NKX family of homeodo...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human NKX2-3
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	38
Uniprot ID	Q8TAU0
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...