



Product Datasheet

MLXIPL antibody (orb67303)


Description	Rabbit polyclonal to MLXIPL. MLX-interacting protein-like (MLXIPL) is a nuclear protein involved in t...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, IP, WB
Immunogen	C-terminal region of mouse MLXIPL
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	95
Uniprot ID	Q99MZ3
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.



Western blot (WB) of total cell extracts...



c) human HeLa, (b



d) human Jurkat, (c) human HepG2; using

...