


Product Datasheet

Histone H2A antibody (orb673037)

Description	Histone H2A antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthesized peptide derived from Human Histone H2A around the non-phosphorylation site of T121.
Target	HIST1H2AG/HIST1H2AB/HIST1H2AD/HIST1H2AH/HIST2H2AA3/HIST3H2A
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04908, Q7L7L0, P0C0S8, Q96KK5, P20671, Q6FI13
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela+UV 5' cell...