


Product Datasheet

PIK3R1 antibody (orb673071)

Description	PIK3R1 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	PIK3R1
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P27986
Dilution Range	WB: 1: 1000-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) 3T3