

Product Datasheet

OR51S1 antibody (orb673220)

Description	OR51S1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 51S1.
Target	OR51S1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NGJ8
Dilution Range	IF: 1:200-1:1000, ELISA: 1:10000
Expiration Date	12 months from date of receipt.