

Product Datasheet

LDLRAD3 antibody (orb673312)

Description	LDLRAD3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Synthesized peptide derived from the Internal region of Human LDLRAD3.
Target	LDLRAD3
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86YD5
Dilution Range	IF: 1:200-1:1000, ELISA: 1:40000
Expiration Date	12 months from date of receipt.