

Product Datasheet

GPBAR1 antibody (orb673325)

Description	GPBAR1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthesized peptide derived from the N-terminal region of Human TGR5.
Target	GPBAR1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TDU6
Dilution Range	WB: 1:500-2000 ELISA: 2000-20000 IHC:1:50-300
Expiration Date	12 months from date of receipt.