

Product Datasheet

GRIN2D antibody (orb673670)

Description	GRIN2D antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from the Internal region of Human NMDAε4.
Target	GRIN2D
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15399
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:40000
Expiration Date	12 months from date of receipt.