

Product Datasheet

OR4K17 antibody (orb673838)

Description	OR4K17 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 4K17.
Target	OR4K17
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NGC6
Dilution Range	IF: 1:200-1:1000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.