


Product Datasheet

GRP78 Antibody (orb67384)

Description	Mouse Anti-Human GRP78 Monoclonal IgG1
Species/Host	Mouse
Reactivity	Bovine, Frog, Fungi, Hamster, Human, Monkey, Mouse, Plant, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	His-tagged human GRP78
Target	GRP78
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	0.5 µg/ml of SMC-196 was sufficient for detection of Grp78 in 10 µg of HeLa cell lysate by ECL immunoblot analysis.
Isotype	IgG1
Clonality	Monoclonal
MW	78kDa
Uniprot ID	P20029
NCBI	NP_001156906.1
Entrez	14828
Dilution Range	WB (1:2000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.



Western blot analysis of GRP78 antibody