

Product Datasheet

PLCG2 (Phospho-Tyr753) antibody (orb6741)

Description Rabbit polyclonal to PLCG2 (Phospho-Tyr753). Enzymes of the phospholipase C family catalyze the hydr...

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated


Tested Applications IHC-P, WB

Note For research use only

Clonality Polyclonal

Dilution Range WB:1:500-1:2000, IHC-P:1:50-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of HepG2 whole cel...