

Product Datasheet

GABA A Receptor Antibody (orb67427)

Description

Mouse Anti-Mouse GABA-A Receptor Beta 3



Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

AM, ICC, IHC, WB

Applications

Immunohistochemical analysis of formalin...

Immunogen

Fusion protein amino acids 370-433 of mouse GABA-A-R-Beta 3

Target

GABA A Receptor

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage


-20°C

Note

For research use only

Application notes

1 µg/ml of SMC-339 was sufficient for detection of Beta 3 GABA receptor in 10 µg of rat brain lysate by colorimetric immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.



Western blot analysis of rat brain membr...

Isotype

IgG1

Clonality

Monoclonal

MW

55kDa

Uniprot ID

[P63080](#)

NCBI

[NP_032097.1](#)

Entrez

[14402](#)

Dilution Range

WB (1:1000), IHC-P (1:1000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.