

Product Datasheet

PTAFR antibody (orb674280)

Description	PTAFR antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Synthesized peptide derived from the Internal region of Human PAF-R.
Target	PTAFR
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P25105
Dilution Range	IF: 1:200-1:1000, ELISA: 1:10000
Expiration Date	12 months from date of receipt.