

Product Datasheet

FOXD4 antibody (orb674487)

Description	FOXD4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from the Internal region of Human FoxD4/D4L.
Target	FOXD4/FOXD4L4/FOXD4L5/FOXD4L3/FOXD4L6
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5VV16 , Q12950 , Q6VB85 , Q8WXT5 , Q6VB84 , Q3SYB3
Dilution Range	WB: 1:500-1:2000, ELISA: 1:40000
Expiration Date	12 months from date of receipt.