

Product Datasheet

BNIP2 antibody (orb674503)

Description	BNIP2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Immunogen	Synthesized peptide derived from the C-terminal region of Human BNIP-2.
Target	BNIP2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q12982
Dilution Range	WB: 1:500-1:2000, IF: 1:200-1:1000, FC: 1:200-1:400, ELISA: 1:40000
Expiration Date	12 months from date of receipt.