


Product Datasheet

VACHT Antibody (orb67481)

Description	Mouse Anti-Human VACHT Monoclonal IgG1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Synthetic peptide amino acids 521-532 of human VACHT
Target	VACHT
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	A dilution of 1:50-1:200 of SMC-341 was sufficient for detection of VACHT Transporter in rat brain using immunohistochemistry analysis and goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	56kDa
Uniprot ID	Q16572
NCBI	NP_003046.2
Entrez	6572
Dilution Range	WB (1:1000), IHC-P (1:200), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.



Western blot analysis of rat membrane ti...