


Product Datasheet

Ankyrin B Antibody (orb67488)

Description	Mouse Anti-Human Ankyrin B Monoclonal IgG1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Synthetic peptide amino acids 203-496 of human Ankyrin-B
Target	Ankyrin B
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	1 µg/ml of SMC-400 was sufficient for detection of Ankyrin-B in 20 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	200kDa
Uniprot ID	Q01484
NCBI	NP_001120965.1
Entrez	287
Dilution Range	WB (1:1000), IHC-P (1:200), ICC/IF (1:200)
Expiration Date	12 months from date of receipt.



Western blot analysis of rat brain membr...