



Product Datasheet

CICN3 Antibody (orb67496)

Description	Mouse Anti-Rat Clcn3 Monoclonal IgG1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Synthetic peptide amino acids 98-115 (cytoplasmic N-terminus) of rat Clcn3
Target	CICN3
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	1 µg/ml of SMC-408 was sufficient for detection of Clcn3 in 20 µg of rat brain membrane lysate and assayed by colorimetric immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	90kDa
Uniprot ID	P51792
NCBI	NP_445815.2
Entrez	84360
Dilution Range	WB (1:1000), IHC-P (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of neuroblas...



Western Blot analysis of rat brain membr...