



## Product Datasheet

### Cleaved-SPTAN1 (D1185) antibody (orb675230)

|                            |                                                                                     |
|----------------------------|-------------------------------------------------------------------------------------|
| <b>Description</b>         | Cleaved-SPTAN1 (D1185) antibody                                                     |
| <b>Species/Host</b>        | Rabbit                                                                              |
| <b>Reactivity</b>          | Human                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                        |
| <b>Tested Applications</b> | ELISA, WB                                                                           |
| <b>Immunogen</b>           | Synthesized peptide derived from the Internal region of Human Spectrin $\alpha$ II. |
| <b>Target</b>              | SPTAN1                                                                              |
| <b>Preservatives</b>       | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.        |
| <b>Form/Appearance</b>     | Liquid                                                                              |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                       |
| <b>Note</b>                | For research use only                                                               |
| <b>Isotype</b>             | IgG                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                          |
| <b>Uniprot ID</b>          | <b>Q13813</b>                                                                       |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, ELISA: 1:10000                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                     |



Western blot analysis of using Cleaved-S...



Western blot analysis of HeLa cells usin...