



## Product Datasheet

### API5 antibody (orb675323)

|                            |                                                                              |
|----------------------------|------------------------------------------------------------------------------|
| <b>Description</b>         | API5 antibody                                                                |
| <b>Species/Host</b>        | Rabbit                                                                       |
| <b>Reactivity</b>          | Human, Mouse                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                 |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                               |
| <b>Immunogen</b>           | Synthesized peptide derived from the C-terminal region of Human API5.        |
| <b>Target</b>              | API5                                                                         |
| <b>Preservatives</b>       | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide. |
| <b>Form/Appearance</b>     | Liquid                                                                       |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                |
| <b>Note</b>                | For research use only                                                        |
| <b>Isotype</b>             | IgG                                                                          |
| <b>Clonality</b>           | Polyclonal                                                                   |
| <b>Uniprot ID</b>          | <a href="#">Q9BZZ5</a>                                                       |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IHC-P: 1:100-1:300, ELISA: 1:20000                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                              |



Western blot analysis of NIH-3T3 cells u...



Western blot analysis of HuvEc cells usi...