

Product Datasheet

Methylated Lysine Antibody (orb67534)

Description	Rabbit Anti-Methylated Lysine Polyclonal
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IP, WB
Immunogen	Methylated KLH Conjugated
Target	Methylated Lysine
Preservatives	PBS, 50% glycerol *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	0.2-0.5 µg/ml of SPC-158 was sufficient for detection of the methylated histone by western blot analysis using melanoma cells in TBSt.
Clonality	Polyclonal
Dilution Range	WB (1:20000), ICC/IF (1:200)
Expiration Date	12 months from date of receipt.



Western blot analysis of Methylated Lysi...