

Product Datasheet

TNFSF9 antibody (orb675395)

Description	TNFSF9 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Synthesized peptide derived from the Internal region of Human CD137L.
Target	TNFSF9
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P41273
Dilution Range	IF: 1:200-1:1000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.