



Product Datasheet

CDYL2 antibody (orb675409)

Description	CDYL2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Synthesized peptide derived from the N-terminal region of Human CDYL2.
Target	CDYL2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N8U2
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:100-1:300, IHC-P: 1:200-1:1000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 Jurkat cells.



Western blot analysis of HEPG2 cells.