



Product Datasheet

TLR4 Antibody (orb67563)

Description	Rabbit Anti-Human TLR4 Polyclonal
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FACS, FC, ICC, IHC, WB
Immunogen	Developed against a synthetic peptide corresponding to amino acids 420-435 of human TLR4
Target	TLR4
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	2 µg/ml of SPC-200 was sufficient for detection of TLR4 in 100 ng of partial recombinant mouse TLR4 protein by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	75-80kDa
Uniprot ID	O00206
NCBI	NP_612564.1
Entrez	7099
Dilution Range	WB (1:500), IHC-P (1:50)
Expiration Date	12 months from date of receipt.



Western blot analysis of TLR4 antibody



Immunohistochemical analysis of formalin...