



Product Datasheet

ENaC Antibody (orb67568)

Description	Rabbit Anti-Rat ENaC beta Polyclonal
Species/Host	Rabbit
Reactivity	Frog, Hamster, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	Produced against the C-terminal tail (amino acids 617-638) of rat beta ENaC (antibody designation 3755-2)
Target	ENaC
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	1 µg/ml of SPC-404 was sufficient for detection of beta-ENaC in 20 µg of rat kidney tissue lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	87kDa
Uniprot ID	P37090
NCBI	NP_036780
Entrez	24767
Dilution Range	WB (1:1000), IHC-P (1:100)
Expiration Date	12 months from date of receipt.



Confocal
Immunofluorescence
analysis ...



Western blot analysis of
rat kidney tiss...