



Product Datasheet

Amyloid Antibody (orb67575)

Description	Rabbit Anti-Human Amyloid Oligomers (A11) Polyclonal
Species/Host	Rabbit
Reactivity	Human, Mouse, Other, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, IP, WB
Immunogen	Synthetic molecular mimic of soluble oligomers
Target	Amyloid Oligomers (A11)
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Storage	-20°C
Note	For research use only
Application notes	A 1:200 dilution of SPC-506 was sufficient for detection of amyloid oligomers in 10 µg of mouse brain lysates by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
Uniprot ID	P05067
NCBI	NM_000484.2
Entrez	8666
Dilution Range	The suggested dilution for the current lot of SPC-506 in Western Blot and Dot Blot is 1:100; optimal dilutions for assays should be determined by the user.
Expiration Date	12 months from date of receipt.



Western blot analysis of Abeta42 fibrils...



Dot blot analysis of Abeta42 and polyQ36...