


Product Datasheet

TUBA1A antibody (orb675946)

Description	TUBA1A antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from Human Tubulin α around the non-acetylation site of K163.
Target	TUBA1A/TUBA1B/TUBA1C/TUBA3C/TUBA4A
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13748 , Q9BQE3 , P68363 , Q71U36 , P68366
Dilution Range	WB: 1:500-1:2000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of K562, PC12 cell...