



Product Datasheet

PRKAA1 antibody (orb676118)

Description	PRKAA1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant Protein
Target	PRKAA1
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13131
Dilution Range	WB: 1: 500-1000, IHC-P: 1: 200-500
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) HeLa, 2) 3T3...



Immunohistochemical staining of Mouse Sk...