


Product Datasheet

Cleaved-PARP1 (D214) antibody (orb676119)


Description	Cleaved-PARP1 (D214) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from the Internal region of Human PARP-1.
Target	PARP1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09874
Dilution Range	WB: 1:500-1:2000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of SH-SY5Y cell lysate.



Western blot analysis of A549 cell lysate.



Western blot analysis of MCF7 and 293T cell lysates.

...